



## SL9C1 Polyclonal Antibody

Cat No: HR1AP11394

For research use only

### Overview

Product Name	SL9C1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	129kDa
GeneID?Human?	SLC9C1 SLC9A10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 9 member C1(SLC9C1) Homo sapiens SLC9A10 is a member of the sodium-hydrogen exchanger (NHE) family (see SLC9A1, MIM 107310) and is required for male fertility and sperm motility (Wang et al., 2003 [PubMed 14634667]).[supplied by OMIM, Apr 2009],