

## SL9C2 Polyclonal Antibody

Cat No: HR1AP11395

For research use only

### Overview

Product Name	SL9C2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	123kDa
GeneID?Human?	SLC9C2 SLC9A11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	domain:Contains an ion transport-like region is related to the membrane segments of voltage-gated ion channels.,function:Involved in pH regulation.,similarity:Belongs to the monovalent cation:proton antiporter 1 (CPA1) transporter (TC 2.A.36) family.,similarity:Contains 1 cyclic nucleotide-binding domain.,