

## SL9A1 Polyclonal Antibody

Cat No: HR1AP11393

For research use only

### Overview

Product Name	SL9A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	89kDa
GeneID?Human?	SLC9A1 APNH1 NHE1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 9 member A1(SLC9A1) Homo sapiens This gene encodes a Na<sup>+</sup>/H<sup>+</sup> antiporter that is a member of the solute carrier family 9. The encoded protein is a plasma membrane transporter that is expressed in the kidney and intestine. This protein plays a central role in regulating pH homeostasis, cell migration and cell volume. This protein may also be involved in tumor growth. [provided by RefSeq, Sep 2011],</p>