

NAC2 Polyclonal Antibody

Cat No: HR1AP11391

For research use only

Overview

Product Name	NAC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	101kDa
GeneID?Human?	SLC8A2 KIAA1087 NCX2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	enzyme regulation:By ATP.,function:Rapidly transports Ca(2+) during excitation-contraction coupling. Ca(2+) is extruded from the cell during relaxation so as to prevent overloading of intracellular stores.,similarity:Belongs to the sodium/potassium/calcium exchanger family. SLC8 subfamily.,similarity:Contains 2 Calx-beta domains.,