

SCNAA Polyclonal Antibody

Cat No: HR1AP11387

For research use only

Overview

Product Name	SCNAA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	215kDa
GeneID?Human?	SCN10A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	sodium voltage-gated channel alpha subunit 10(SCN10A) Homo sapiens The protein encoded by this gene is a tetrodotoxin-resistant voltage-gated sodium channel alpha subunit. The properties of the channel formed by the encoded transmembrane protein can be altered by interaction with different beta subunits. This protein may be involved in the onset of pain associated with peripheral neuropathy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014],