

## S6A18 Polyclonal Antibody

Cat No: HR1AP11383

For research use only

### Overview

Product Name	S6A18 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	69kDa
GeneID?Human?	SLC6A18 XTRP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 6 member 18(SLC6A18) Homo sapiens The SLC6 family of proteins, which includes SLC6A18, act as specific transporters for neurotransmitters, amino acids, and osmolytes like betaine, taurine, and creatine. SLC6 proteins are sodium cotransporters that derive the energy for solute transport from the electrochemical gradient for sodium ions (Hoglund et al., 2005 [PubMed 16125675]).[supplied by OMIM, Apr 2010],</p>