

SCNM1 Polyclonal Antibody

Cat No: HR1AP11386

For research use only

Overview

Product Name	SCNM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	SCNM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	sodium channel modifier 1(SCNM1) Homo sapiens SCNM1 is a zinc finger protein and putative splicing factor. In mice, Scnm1 modifies phenotypic expression of Scn8a (MIM 600702) mutations (Buchner et al., 2003 [PubMed 12920299]).[supplied by OMIM, Oct 2009],