

S6A13 Polyclonal Antibody

Cat No: HR1AP11382

For research use only

Overview

Product Name	S6A13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	SLC6A13 GAT2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Terminates the action of GABA by its high affinity sodium-dependent reuptake into presynaptic terminals.,similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,tissue specificity:Expressed in brain, kidney, lung, liver and testis.,