

## S6A12 Polyclonal Antibody

Cat No: HR1AP11381

For research use only

### Overview

Product Name	S6A12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	SLC6A12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Transports betaine and GABA. May have a role in regulation of GABAergic transmission in the brain through the reuptake of GABA into presynaptic terminals, as well as in osmotic regulation.,similarity:Belongs to the sodium:neurotransmitter symporter (SNF) family.,subunit:Interacts with LIN7C.,tissue specificity:Liver, heart, skeletal muscle, placenta, and a widespread distribution in the brain.,