

# DGLA Polyclonal Antibody

Cat No: HR1AP11380

For research use only

## Overview

Product Name	DGLA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	114kDa
GeneID?Human?	DAGLA C11orf11 KIAA0659 NSDDR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	diacylglycerol lipase alpha(DAGLA) Homo sapiens This gene encodes a diacylglycerol lipase. The encoded enzyme is involved in the biosynthesis of the endocannabinoid 2-arachidonoyl-glycerol.[provided by RefSeq, Nov 2010],