



## SSA27 Polyclonal Antibody

Cat No: HR1AP11372

For research use only

### Overview

Product Name	SSA27 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	SSSCA1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Sjogren syndrome/scleroderma autoantigen 1(SSSCA1) Homo sapiens This antigen is recognized by a subset of anti-centromere antibodies from patients with scleroderma and/or Sjogren's syndrome. Subcellular localization has not yet been established. [provided by RefSeq, Jul 2008],