

TRPC1 Polyclonal Antibody

Cat No: HR1AP11368

For research use only

Overview

Product Name	TRPC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	TRPC1 TRP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>transient receptor potential cation channel subfamily C member 1(TRPC1) Homo sapiens The protein encoded by this gene is a membrane protein that can form a non-selective channel permeable to calcium and other cations. The encoded protein appears to be induced to form channels by a receptor tyrosine kinase-activated phosphatidylinositol second messenger system and also by depletion of intracellular calcium stores. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011],</p>