

SASH3 Polyclonal Antibody

Cat No: HR1AP11365

For research use only

Overview

Product Name	SASH3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	SASH3 CXorf9 SLY
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SAM and SH3 domain containing 3(SASH3) Homo sapiens The protein encoded by this gene contains a Src homology-3 (SH3) domain and a sterile alpha motif (SAM), both of which are found in proteins involved in cell signaling. This protein may function as a signaling adapter protein in lymphocytes.[provided by RefSeq, Sep 2009],