

## SH3K1 Polyclonal Antibody

Cat No: HR1AP11364

For research use only

### Overview

Product Name	SH3K1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	SH3KBP1 CIN85
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SH3 domain containing kinase binding protein 1(SH3KBP1) Homo sapiens This gene encodes an adapter protein that contains three N-terminal Src homology domains, a proline rich region and a C-terminal coiled-coil domain. The encoded protein facilitates protein-protein interactions and has been implicated in numerous cellular processes including apoptosis, cytoskeletal rearrangement, cell adhesion and in the regulation of clathrin-dependent endocytosis. Alternate splicing results in multiple transcript variants.[provided by RefSeq, May 2010],