

SH3B4 Polyclonal Antibody

Cat No: HR1AP11363

For research use only

Overview

Product Name	SH3B4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	105kDa
GeneID?Human?	SH3BP4 BOG25 EHB10 TTP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SH3 domain binding protein 4(SH3BP4) Homo sapiens This gene encodes a protein with 3 Asn-Pro-Phe (NPF) motifs, an SH3 domain, a PXXP motif, a bipartite nuclear targeting signal, and a tyrosine phosphorylation site. This protein is involved in cargo-specific control of clathrin-mediated endocytosis, specifically controlling the internalization of a specific protein receptor. [provided by RefSeq, Jul 2008],