

3BP1 Polyclonal Antibody

Cat No: HR1AP11361

For research use only

Overview

Product Name	3BP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	SH3BP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Binds differentially to the SH3 domains of certain proteins of signal transduction pathways. This protein binds preferentially to the ABL1 proto-oncogene, SRC and GRB2. Shows GAP activity for Rac-related proteins but not for Rho- or Ras-related proteins. It inhibits PDGF-induced membrane ruffling mediated by Rac.,sequence caution:Translated as Cys.,similarity:Contains 1 BAR domain.,similarity:Contains 1 Rho-GAP domain.,