

SH3L3 Polyclonal Antibody

Cat No: HR1AP11359

For research use only

Overview

Product Name	SH3L3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	SH3BGRL3 P1725
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Could act as a modulator of glutaredoxin biological activity.,similarity:Belongs to the SH3BGR family.,similarity:Contains 1 glutaredoxin domain.,tissue specificity:Ubiquitous.,