

SH3L2 Polyclonal Antibody

Cat No: HR1AP11358

For research use only

Overview

Product Name	SH3L2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	SH3BGRL2 FASH3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the SH3BGR family.,tissue specificity:Highly expressed in brain, placenta, liver and kidney. Expressed in retina.,