

STAC3 Polyclonal Antibody

Cat No: HR1AP11356

For research use only

Overview

Product Name	STAC3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	STAC3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SH3 and cysteine rich domain 3(STAC3) Homo sapiens The protein encoded by this gene is a component of the excitation-contraction coupling machinery of muscles. This protein is a member of the Stac gene family and contains an N-terminal cysteine-rich domain and two SH3 domains. Mutations in this gene are a cause of Native American myopathy. [provided by RefSeq, Nov 2013],