



STAC Polyclonal Antibody

Cat No: HR1AP11355

For research use only

Overview

Product Name	STAC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	STAC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Probably involved in a neuron-specific signal transduction.,similarity:Contains 1 phorbol-ester/DAG-type zinc finger.,similarity:Contains 1 SH3 domain.,