

SH23A Polyclonal Antibody

Cat No: HR1AP11351

For research use only

Overview

Product Name	SH23A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GeneID?Human?	SH2D3A NSP1 UNQ175/PRO201
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in JNK activation.,PTM:Phosphorylated on tyrosine.,similarity:Contains 1 SH2 domain.,subunit:Interacts with BCAR1.,tissue specificity:Weakly expressed in placenta, fetal kidney, fetal lung, adult pancreas, adult kidney and adult lung.,