

SCAFB Polyclonal Antibody

Cat No: HR1AP11348

For research use only

Overview

Product Name	SCAFB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	160kDa
GeneID?Human?	SCAF11 CASP11 SFRS2IP SIP1 SRSF2IP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Plays a role in pre-mRNA alternative splicing by regulating spliceosome assembly.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Interacts with SFRS2, U2AF2 and SNRNP70.,tissue specificity:Widely expressed.,