



SPI2 Polyclonal Antibody

Cat No: HR1AP11347

For research use only

Overview

Product Name	SPI2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	SERPINI2 MEPI PI14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>serpin family I member 2(SERPINI2) Homo sapiens The gene encodes a member of a family of proteins that acts as inhibitors of serine proteases. These proteins function in the regulation of a variety of physiological processes, including coagulation, fibrinolysis, development, malignancy, and inflammation. Expression of the encoded protein may be downregulated during pancreatic carcinogenesis. Alternative splicing results in multiple transcript variants for this gene. [provided by RefSeq, Jan 2013],</p>