

## SPA9 Polyclonal Antibody

Cat No: HR1AP11344

For research use only

### Overview

Product Name	SPA9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	SERPINA9 GCET1 SERPINA11 UNQ692/PRO1337
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the serpin family.,tissue specificity:Highly expressed in normal germinal center (GC) B-cells and GC B-cell-derived malignancies.,