

CNEP1 Polyclonal Antibody

Cat No: HR1AP11341

For research use only

Overview

Product Name	CNEP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	CTDNEP1 DULLARD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:A phosphoprotein + H(2)O = a protein + phosphate.,function:Serine/threonine phosphatase which may be required for proper nuclear membrane morphology. Involved in LPIN1 dephosphorylation. May antagonize BMP signaling.,similarity:Belongs to the dullard family.,similarity:Contains 1 FCP1 homology domain.,