

PP4R1 Polyclonal Antibody

Cat No: HR1AP11340

For research use only

Overview

Product Name	PP4R1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	104kDa
GeneID?Human?	PPP4R1 MEG1 PP4R1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>protein phosphatase 4 regulatory subunit 1(PPP4R1) Homo sapiens This gene encodes one of several alternate regulatory subunits of serine/threonine protein phosphatase 4 (PP4). The protein features multiple HEAT repeats. This protein forms a complex with PP4RC. This complex may have a distinct role from other PP4 complexes, including regulation of HDAC3 (Zhang et al., PMID: 15805470). There is also a transcribed pseudogene on chromosome 20. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2012],</p>