

## PP2BC Polyclonal Antibody

Cat No: HR1AP11339

For research use only

### Overview

Product Name	PP2BC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	PPP3CC CALNA3 CNA3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>protein phosphatase 3 catalytic subunit gamma(PPP3CC) Homo sapiens Calcineurin is a calcium-dependent, calmodulin-stimulated protein phosphatase involved in the downstream regulation of dopaminergic signal transduction. Calcineurin is composed of a regulatory subunit and a catalytic subunit. The protein encoded by this gene represents one of the regulatory subunits that has been found for calcineurin. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],</p>