

## PP1RA Polyclonal Antibody

Cat No: HR1AP11332

For research use only

### Overview

Product Name	PP1RA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	103kDa
GeneID?Human?	PPP1R10 CAT53 FB19 PNUTS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>protein phosphatase 1 regulatory subunit 10(PPP1R10) Homo sapiens This gene encodes a protein phosphatase 1 binding protein. The encoded protein plays a role in many cellular processes including cell cycle progression, DNA repair and apoptosis by regulating the activity of protein phosphatase 1. This gene lies within the major histocompatibility complex class I region on chromosome 6, and alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012],</p>