

## SGK2 Polyclonal Antibody

Cat No: HR1AP11326

For research use only

### Overview

Product Name	SGK2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	SGK2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SGK2, serine/threonine kinase 2(SGK2) Homo sapiens This gene encodes a serine/threonine protein kinase. Although this gene product is similar to serum- and glucocorticoid-induced protein kinase (SGK), this gene is not induced by serum or glucocorticoids. This gene is induced in response to signals that activate phosphatidylinositol 3-kinase, which is also true for SGK. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2010],</p>