

PRP4B Polyclonal Antibody

Cat No: HR1AP11323

For research use only

Overview

Product Name	PRP4B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	110kDa
GeneID?Human?	PRPF4B KIAA0536 PRP4 PRP4H PRP4K
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	pre-mRNA processing factor 4B(PRPF4B) Homo sapiens Pre-mRNA splicing occurs in two sequential transesterification steps, and the protein encoded by this gene is thought to be involved in pre-mRNA splicing and in signal transduction. This protein belongs to a kinase family that includes serine/arginine-rich protein-specific kinases and cyclin-dependent kinases (CDKs). This protein is regarded as a CDK-like kinase (Clk) with homology to mitogen-activated protein kinases (MAPKs). [provided by RefSeq, Jul 2008],