

## PLK4 Polyclonal Antibody

Cat No: HR1AP11322

For research use only

### Overview

Product Name	PLK4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	106kDa
GeneID?Human?	PLK4 SAK STK18
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>polo like kinase 4(PLK4) Homo sapiens This gene encodes a member of the polo family of serine/threonine protein kinases. The protein localizes to centrioles, complex microtubule-based structures found in centrosomes, and regulates centriole duplication during the cell cycle. Three alternatively spliced transcript variants that encode different protein isoforms have been found for this gene. [provided by RefSeq, Jun 2010],</p>