

NEK1 Polyclonal Antibody

Cat No: HR1AP11318

For research use only

Overview

Product Name	NEK1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	138kDa
GeneID?Human?	NEK1 KIAA1901
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>NIMA related kinase 1(NEK1) Homo sapiens The protein encoded by this gene is a serine/threonine kinase involved in cell cycle regulation. The encoded protein is found in a centrosomal complex with FEZ1, a neuronal protein that plays a role in axonal development. Defects in this gene are a cause of polycystic kidney disease (PKD). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2010],</p>