

NEK10 Polyclonal Antibody

Cat No: HR1AP11319

For research use only

Overview

Product Name	NEK10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	128kDa
GeneID?Human?	NEK10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,cofactor:Magnesium.,similarity:Belongs to the protein kinase superfamily. NEK Ser/Thr protein kinase family. NIMA subfamily.,similarity:Contains 1 ARM repeat.,similarity:Contains 1 protein kinase domain.,