

STK31 Polyclonal Antibody

Cat No: HR1AP11313

For research use only

Overview

Product Name	STK31 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	112kDa
GeneID?Human?	STK31 SGK396
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serine/threonine kinase 31(STK31) Homo sapiens This gene is similar to a mouse gene that encodes a putative protein kinase with a tudor domain, and shows testis-specific expression. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],