



## PRS23 Polyclonal Antibody

Cat No: HR1AP11309

For research use only

### Overview

Product Name	PRS23 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	PRSS23 ZSIG13 UNQ270/PRO307
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	protease, serine 23(PRSS23) Homo sapiens This gene encodes a conserved member of the trypsin family of serine proteases. Mouse studies found a decrease of mRNA levels of this gene after ovulation was induced. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014],