



SFRP3 Polyclonal Antibody

Cat No: HR1AP11306

For research use only

Overview

Product Name	SFRP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	FRZB FIZ FRE FRP FRZB1 SFRP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	frizzled-related protein(FRZB) Homo sapiens The protein encoded by this gene is a secreted protein that is involved in the regulation of bone development. Defects in this gene are a cause of female-specific osteoarthritis (OA) susceptibility. [provided by RefSeq, Apr 2010],