

SCAR3 Polyclonal Antibody

Cat No: HR1AP11304

For research use only

Overview

Product Name	SCAR3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	SCARA3 CSR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	scavenger receptor class A member 3(SCARA3) Homo sapiens This gene encodes a macrophage scavenger receptor-like protein. This protein has been shown to deplete reactive oxygen species, and thus play an important role in protecting cells from oxidative stress. The expression of this gene is induced by oxidative stress. Alternatively spliced transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Jul 2008],