

# RPAP3 Polyclonal Antibody

Cat No: HR1AP11300

For research use only

## Overview

Product Name	RPAP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	RPAP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RNA polymerase II associated protein 3(RPAP3) Homo sapiens This gene encodes an RNA polymerase II-associated protein. The encoded protein may function in transcriptional regulation and may also regulate apoptosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009],