

RPAP1 Polyclonal Antibody

Cat No: HR1AP11299

For research use only

Overview

Product Name	RPAP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	153kDa
GeneID?Human?	RPAP1 KIAA1403
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RNA polymerase II associated protein 1(RPAP1) Homo sapiens This protein forms part of the RNA polymerase II (RNAPII) enzyme complex and may recruit RNAPII to chromatin through its interaction with acetylated histones. [provided by RefSeq, Jul 2012],