

RPF1 Polyclonal Antibody

Cat No: HR1AP11295

For research use only

Overview

Product Name	RPF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	RPF1 BXDC5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be required for ribosome biogenesis.,similarity:Contains 1 Brix domain.,