

BMS1 Polyclonal Antibody

Cat No: HR1AP11293

For research use only

Overview

Product Name	BMS1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	141kDa
GeneID?Human?	BMS1 BMS1L KIAA0187
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	BMS1, ribosome biogenesis factor(BMS1) Homo sapiens This gene likely encodes a ribosome assembly protein. A similar protein in yeast functions in 35S-rRNA processing, which includes a series of cleavage steps critical for formation of 40S ribosomes. Related pseudogenes exist on chromosomes 2, 9, 10, 15, 16, and 22.[provided by RefSeq, Mar 2009],