

RHXF2 Polyclonal Antibody

Cat No: HR1AP11290

For research use only

Overview

Product Name	RHXF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	RHOXF2 PEPP2 THG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Rhox homeobox family member 2(RHOXF2) Homo sapiens This gene, which encodes a transcriptional repressor, is one of two paralogous X-linked homeobox-containing genes and is highly expressed in a variety of cancers. In addition, the encoded protein associates with the cell membrane and with microtubules, and is concentrated at the leading edge of migratory cells. [provided by RefSeq, Dec 2015],