

# RHXF1 Polyclonal Antibody

Cat No: HR1AP11289

For research use only

## Overview

Product Name	RHXF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	RHOXF1 OTEX PEPP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Rhox homeobox family member 1(RHOXF1) Homo sapiens This gene is a member of the PEPP subfamily of paired-like homeobox genes. The gene may be regulated by androgens and epigenetic mechanisms. The encoded nuclear protein is likely a transcription factor that may play a role in human reproduction. [provided by RefSeq, Dec 2012],