



RHOF Polyclonal Antibody

Cat No: HR1AP11285

For research use only

Overview

Product Name	RHOF Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RHOF ARHF RIF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Plasma membrane-associated small GTPase which cycles between an active GTP-bound and an inactive GDP-bound state. Causes the formation of thin, actin-rich surface projections called filopodia. Functions cooperatively with CDC42 and Rac to generate additional structures, increasing the diversity of actin-based morphology.,similarity:Belongs to the small GTPase superfamily. Rho family.,