



RHOJ Polyclonal Antibody

Cat No: HR1AP11286

For research use only

Overview

Product Name	RHOJ Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RHOJ ARHJ RASL7B RHOI TCL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ras homolog family member J(RHOJ) Homo sapiens This gene encodes one of the many small GTP-binding proteins in the Rho family shown to be associated with focal adhesions in endothelial cells (PMID: 21148427, 22103495). The encoded protein is activated by vascular endothelial growth factor and may regulate angiogenesis. [provided by RefSeq, Dec 2011],