

RND3 Polyclonal Antibody

Cat No: HR1AP11284

For research use only

Overview

Product Name	RND3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	RND3 ARHE RHO8 RHOE
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Rho family GTPase 3(RND3) Homo sapiens This gene encodes a protein which is a member of the small GTPase protein superfamily. The encoded protein binds only GTP but has no GTPase activity, and appears to act as a negative regulator of cytoskeletal organization leading to loss of adhesion. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Dec 2011],