

RND1 Polyclonal Antibody

Cat No: HR1AP11283

For research use only

Overview

Product Name	RND1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	RND1 RHO6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Rho family GTPase 1(RND1) Homo sapiens This gene encodes a protein that belongs to the Rho GTPase family. Members of this family regulate the organization of the actin cytoskeleton in response to extracellular growth factors. A similar protein in rat interacts with a microtubule regulator to control axon extension. [provided by RefSeq, Apr 2014],