

# RHPN2 Polyclonal Antibody

Cat No: HR1AP11282

For research use only

## Overview

Product Name	RHPN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	RHPN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>rhophilin Rho GTPase binding protein 2(RHPN2) Homo sapiens This gene encodes a member of the rhophilin family of Ras-homologous (Rho)-GTPase binding proteins. The encoded protein binds both GTP- and GDP-bound RhoA and GTP-bound RhoB and may be involved in the organization of the actin cytoskeleton. [provided by RefSeq, Apr 2009],</p>