

ARHG4 Polyclonal Antibody

Cat No: HR1AP11280

For research use only

Overview

Product Name	ARHG4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	ARHGEF4 KIAA1112
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Rho guanine nucleotide exchange factor 4 (ARHGEF4) Homo sapiens Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein coupled receptors. The protein encoded by this gene may form complex with G proteins and stimulate Rho-dependent signals. Multiple alternatively spliced transcript variants encoding different isoforms have been found, but the full-length nature of some variants has not been determined. [provided by RefSeq, Jun 2013].</p>