

RHG44 Polyclonal Antibody

Cat No: HR1AP11278

For research use only

Overview

Product Name	RHG44 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	89kDa
GeneID?Human?	ARHGAP44 KIAA0672 RICH2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state. Acts as a GTPase activator in vitro for CDC42 and RAC1.,sequence caution:Translation N-terminally extended.,similarity:Contains 1 BAR domain.,similarity:Contains 1 Rho-GAP domain.,tissue specificity:Highly expressed in brain. Expressed at weak level in other tissues.,